Computerized Model Fitting Approach for the NMR Analysis of Polymers[†]

H. N. CHENG

Research Center, Hercules Incorporated, Wilmington, Delaware 19894

Received May 21, 1986; Revised Manuscript Received October 29, 1986

A computer program (called FITCO) was developed for the general analysis of the NMR spectra of polymers through the (analytical) model-fitting approach. The method attempts to optimize the amount of information available from each spectrum and is especially useful for spectra with overlapping resonances and for studies of polymerization mechanisms. Examples are shown of the use of the program for the analysis of copolymer composition in ethyl acrylate-methyl methacrylate copolymers and for the studies of tacticity in polypropylene.

INTRODUCTION

NMR spectroscopy has been widely used for studies of polymer microstructure. Information available includes homopolymer tacticity, copolymer sequence determination, and polymer chain branching. Over the years numerous papers and texts¹⁻⁵ have appeared on this subject.

The analysis of polymer NMR spectra usually follows familiar patterns. A schematic is given in Figure 1 showing the logical steps needed. In general, NMR analysis of the polymer spectrum can be carried out at several levels. At the simplest level, one can treat the polymer spectrum as a fingerprint pattern and use it, for example, for the identification of unknown samples. If more information is needed about the polymer, then the spectrum must be interpreted. This process may be aided by model components or analogous polymers. Recently many advanced techniques⁶⁻⁹ (e.g., polarization transfer, spectral editing, and two-dimensional NMR) have been developed and increasingly used in polymer NMR.9 Assuming that all the significant resonances can be properly assigned, one must then devise computational schemes to obtain polymer compositions and sequence distributions. This general method of analysis can be referred to as the "analytical approach". (An alternative method, called the "synthetic approach", has also been developed 10,11 recently).

A more refined analytical approach is to approximate the copolymerization reaction with a statistical model (Figure 1).¹² One can then associate every spectral intensity with a theoretical expression involving reaction probability parameters. The observed and the theoretical intensities for all the spectral lines are then compared, and optimization is carried out to obtain the best-fit values of the reaction probability parameters. Depending on the goodness of fit, these probability parameters then fully describe the structure of the polymer system in question. This model-fitting approach has been successfully applied to a number of specific polymer systems. ^{12–17,29}

In this work a general computer program has been developed that is capable of rapidly and conveniently applying this model-fitting approach to a variety of polymer systems. Two examples are shown that illustrate this approach: compositional analysis of ethyl acrylate-methyl methacrylate copolymers and determination of polypropylene tacticity by ¹³C NMR. This is the first time that the ¹³C NMR spectra of ethyl acrylate-methyl methacrylate copolymers have been assigned and analyzed.

PROGRAM FITCO

The program is organized into eight sections. The relationships among the various sections are summarized in Figure 2. These sections are described as follows:

[†]Presented at the 190th National Meeting of the American Chemical Society, Chicago, IL, September 1985 (Paper CINF 26).

- (1) The first section involves input of observed spectral intensities. It also produces a header in the output.
- (2) In the second section, the user may wish to devise computational schemes to calculate polymer composition and sequence distribution.
- (3) The format of the input data (in section 1) is specified in this section. One method is to use DATA-READ statements; corrections on the input can be easily made (e.g., lines 200-420, Figure 3).
- (4) The expressions for reaction probability parameters are coded in this section. The common statistical models are (a) Bernoullian, (b) first-order Markov, (c) second-order Markov, (d) Coleman–Fox, and (e) enantiomorphic models. Other models have also been proposed. For the common statistical models the theoretical expressions for many stereosequences have been given in the literature. 18-20 Expressions that are not available may be derived if needed.
- (5) The simplex algorithm is contained here. The purpose is to compare the observed vs. the calculated intensities and to provide a fast and logical means to obtain the optimal values of the reaction probability parameters.
- (6) This section provides the output, viz., final reaction probabilities, optimal values for polymer composition and sequence distribution, and reactivity ratios product.
- (7) This section provides additional opportunities for computations.
- (8) Termination occurs last with provision for looping back to carry out other computations.

It is clear that sections 2, 3, 4, and 7 must be inserted by the user. The other sections need no user input.

EXAMPLES

Program FITCO in its present form is very versatile. It can be used for the compositional analysis and sequence distribution of copolymers and also for the tacticity of homopolymers. Examples are given here to illustrate these cases.

1. Compositional Analysis. The ¹³C NMR spectrum of a copolymer of ethyl acrylate (EA) and methyl methacrylate (MMA) is shown in Figure 4. The complex pattern indicates the combined effects of copolymer sequence placements and tacticity on the ¹³C shifts. In polymers of such complexity, complete interpretation of all resolvable resonances is very difficult and has been accomplished in very few cases. For EA-MMA copolymers, the interpretation of ¹³C NMR spectra has not been previously reported. A close scrutiny of Figure 4, however, indicates that the backbone carbons (α and β) all resonate in the 33-55 ppm region. Except for this region, all other resonances can be assigned with some effort. The assignments were aided by comparing copolymer samples with different compositions. The results are shown in Table I and Figure 4. Of particular interest are the assignments of the methyl region (15-23 ppm) and the carbonyl region (173-180 ppm), where the effects of comonomer sequence placements




Figure 1. Logical steps in the analytical approach for the analysis of NMR spectra of polymers.




Figure 2. Schematic diagram of the FITCO program.

Table I. Assignments of the ¹³C NMR Spectra of Ethyl Acrylate-Methyl Methacrylate Copolymers

designation ^a	shifts, ppm	assignments ^b			
a ₁	60.2	Z (OCH ₂ CH ₃)			
a ₂	51.6	Y (OCH ₃)			
a ₃	13.8	$Z (OCH_2CH_3)$			
b ₁	19.3-23.0	YYY (mm), YYZ (m), ZYZ			
b_2	17.3-19.3	YYY (mr), YYZ (r)			
b_3	15.6-17.3	YYY (rr)			
c _i	179.0-177.4	YYY (rr)			
c_2	177.4-176.0	YYY (mr), YYZ (r),, YZY			
c ₃	176.0-175.0	YYY (mm), YYZ (m),, ZZY			
C ₄	175.0-174.4	ZYZ,			
C ₅	174.4-173.5	, ZZZ			

 $[^]ab$ corresponds to the methyl resonances, and c to the carboxylate resonances. bY = methyl methacrylate; Z = ethyl acrylate.

and tacticity have been unscrambled.

Information on copolymer composition is easy to obtain. One can simply take the ratio of peaks $b_1 + b_2 + b_3$ and a_3 , for example. However, if one needs more information from this partially overlapped spectrum, then the computerized reaction probability model (FITCO) may be used. In this case one may use the Bernoullian statistical model. Let Y = MMA and Z = EA. The Bernoullian probabilities are y_m , y_r , and z, where y_m and y_r are the probabilities of the addition of meso MMA and racemic MMA to a propagating polymer radical, and z is the probability of the addition of EA. The ¹³C NMR spectra of peaks a_1 and a_3 are not sensitive enough to the tacticity of EA to permit separate determination of z_m and z_r .

The theoretical expressions for the Bernoullian probabilities are given in Table II. These theoretical expressions are coded

Table II. Theoretical Bernoullian Expressions and Observed Intensities for Three Samples of EA-MMA Copolymers

designation		obsd intensities ^b			
	theoretical expression ^a	1	2	3	
a ₁	z	20.0	50.5	81.6	
a_2	$y_m + y_r$	80.0	49.5	18.4	
a_3	z	20.2	50.1	78.0	
\mathfrak{b}_1	$k_1(y_m+z)^2$	19.8	24.1	17.9	
b_2	$2k_1y_r(y_m+z)$	40.6	24.9	6.4	
b ₃	$k_1 v_r^2$	27.6	8.4	2.1	
c ₁	$k_2 y_r^2$	25.0	7.4	0.8	
c_2	$2k_{2}v_{m}v_{r} + 2k_{2}v_{r}z + k_{3}^{2}z$	51.6	33.3	6.5	
c ₃	$k_2 v_m^2 + 2k_2 v_m z + 2k_2 z^2$	18.4	35.5	28.1	
C ₄	$k_2 z^2$	3.0	11.0	13.7	
c ₅	$ \begin{array}{l} 2k_2 y_m y_r + 2k_2 y_r z + k_2^2 z \\ k_2 y_m^2 + 2k_2 y_m z + 2k_2 z^2 \\ k_2 z^2 \\ z^3 \end{array} $	1.0	12.7	50.8	

 $^ak_1 = (y_m + y_r)/(y_m + y_r + z)^2$; $k_2 = (y_m + y_r)$; a constant multiplicative factor of 100 is used for all the expressions. b The intensities of the carboxy carbons (c_i) have been separately normalized to 100.

into the program FITCO (lines 3000-3999, Figure 2). The suitable input format statements are incorporated (lines 2000-2999, Figure 2). The program is ready now to accept input. A full listing of this program, written for the EA-MMA copolymer system is shown in Figure 3.

The observed intensities for three samples of EA-MMA copolymers are given in Table II. These values are directly entered into the program FITCO. Initial guess values are provided to start the simplex optimization. The program then takes over and gives the following results:

sample	y_m	y_r	z
1	0.232	0.176	0.053
2	0.578	0.322	0.164
3	0.201	0.502	0.782

Thus in all cases meso/racemic = 0.30/0.70 for the MMA tacticity. This information would be difficult and tedious to obtain by hand calculation. The use of a computerized reaction probability model thus simplifies the analysis and (by using all the observable spectral intensities) decreases the experimental error.

2. Homopolymer Tacticity. The tacticity of polypropylene is a much studied problem. The NMR assignments are well-known at present, 9b,21-23 and attempts have been made to fit the observed tacticities to different statistical models of stereospecific polymerization. At various times, Bernoullian, 22,25 Markovian, 11,22,25 and enantiomorphic site 25,26 models have been used. The latest results seem to indicate that the bicatalytic site model 24,25 is appropriate.

Program FITCO can readily handle this problem and carry out the model fitting. In this case the meso (m) and the racemic (r) additions of the propylene units can be treated like copolymerization. For the first-order Markovian process, two probabilities determine the polymer structure: P_{mr} and P_{rm} . For the second order, four parameters are needed $(\alpha, \beta, \gamma, \delta)$, and for the bicatalytic site model, 25 three parameters (β, σ, ω) . The expressions for Markovian and bicatalytic site probabilities are given in Table III.

The methyl region of the 13 C NMR spectrum has been generally used for tacticity measurements of polypropylene. A spectrum at 90 MHz is given in Figure 5. From the observed methyl pentad intensities, the FITCO program again can readily provide the desired reaction probabilities. The data for several samples taken from the literature are given in Table IV. Sample A resulted from the epimerization studies of Suter and Neuenschwander. The data fitted the Markovian model reasonably well; the second order gave a slightly better fit than the first. As expected, P_{mr} and P_{rm} are both close to 0.5, corresponding to the totally random case. R is the response factor in the simplex algorithm. It is taken to be the mean derivation between the calculated and the observed pentad

```
4835 FOR I = J + 1 TO N9
4840 IF U(J) ( = U(I) THEN 4875
4845 LET D9 = U(I)
4850 LET U(I) = U(J)
4855 LET U(J) = D9
4860 LET D9 = M(I)
4860 LET M(J) = D9
4875 NEXT I
4880 NEXT J
                 REM PROGRAM FITCO/EAMMA
REM DESIGNED FOR THE C-13 NMR ANALYSIS OF
ETHYL ACRYLATE/METHYL METHACRYLATE COPOLYMERS
REM WRITTEN BY H. N. CHENG, HERCULES RESEARCH CENTER,
REM WILMINGTON, DELAWARE 19894, USA
                 REM
                  REM

REM (FOR 11-AREA FIT)

DIM R(9),S(10,8),D(9),Y(8),X(8)

DIM U(8),M(8),L(8),G(40),H(40)

LET V9 = 11

LET T = 0

PRINT "DATE = ";

INPUT A$

PRINT "SAMPLE DESIGNATION = ";

INPUT B$

PRINT "SAMPLE DESIGNATION = ";

PRINT "SAMPLE DESIGNATION = ";

PRINT "SAMPLE DESIGNATION = ";

PRINT "SAMPLE DESIGNATION = ";
     30
                                                                                                                                                                                                                                                                                                                    4880
4890
4900
                                                                                                                                                                                                                                                                                                                  4880 NEXT J
4890 RETURN
4900 REM THIS SUBROUTINE CALCULATES THE CENTROID P9
4905 FOR J = 1 TO NS
4910 LET O(J) = 0
4915 FOR I = 1 TO NS
4916 FOR K = 1 TO NS
4917 IF M(I) = L(K) THEN 4920
4918 NEXT K
4920 LET D(I) = 5(K I) + 0(I)
      130
      160
170
                       PRINT
LET I = 0
READ V1$, V2$, C
LET I = I + 1
IF I ( T THEN 210
PRINT I;"-"; TAB( 6); V1$; V2$; " AREA (";C; " PPM) = ";
INPUT H(I)
IF H(I) ( 0 THEN 400
IF H(I) ) 1000 THEN 480
IF I = V9 THEN 500
GOTO 210
LET T = I + H(I)
IF T ( 0 THEN 480
RESTORE
GOTO 200 .
     215
                                                                                                                                                                                                                                                                                                                    4920
                                                                                                                                                                                                                                                                                                                                         LET D(J) = B(K,J) + D(J)
      250
                                                                                                                                                                                                                                                                                                                    4925
                                                                                                                                                                                                                                                                                                                   4930
4935
4940
                                                                                                                                                                                                                                                                                                                                           LET D(J) = D(J) / N8
     240
250
                                                                                                                                                                                                                                                                                                                                         RETURN
                                                                                                                                                                                                                                                                                                                                         RETURN
REM SUBROUTINE TO PRINT OUT DATA
LET N7 = N7 + 1
PRINT N7; TAB( 4);N6; TAB( 9);
FOR I = 1 TO N8
PRINT M(I);" ";
                                                                                                                                                                                                                                                                                                                    4950
4955
      290
                                                                                                                                                                                                                                                                                                                    4960
      400
                                                                                                                                                                                                                                                                                                                    4965
                                                                                                                                                                                                                                                                                                                   4970
4975
4980
                                                                                                                                                                                                                                                                                                                                        NEXT I
PRINT TAB( 25);10 • R(9); TAB( 37);
      420
                        RESTORE
GOTO 200 .

PRINT "WRONG ENTRY. PLEASE CHECK."
GOTO 230
REM NORMALIZE THE INTENSITIES
LET H9 = H(4) + H(5)
IF H9 = 100 THEN 590
FOR I = 1 TO 6
LET H(I) = H(I) • 100 / H9
NEYT I
      430
                                                                                                                                                                                                                                                                                                                                        FOR J = 1 TO NB

LET S9 = 0.00001 • INT (100000 * S(J9,J))

PRINT S9;" ";
                                                                                                                                                                                                                                                                                                                   4985
4990
      500
                                                                                                                                                                                                                                                                                                                    4991
                                                                                                                                                                                                                                                                                                                   4995
4996
4998
                                                                                                                                                                                                                                                                                                                                         NEXT J
PRINT TAB( 75);"(";Z$;")"
GOSUB 6200
      560
                                                                                                                                                                                                                                                                                                                                         RETURN
REM SIMPLEX SUBROUTINE
PRINT
                                                                                                                                                                                                                                                                                                                    4999
                        LET H(I) = H(I) • 100 / H9

NEXT I

LET H9 = 0

FOR I = 7 TO 11

LET H9 = H9 + H(I)

NEXT I

FOR I = 7 TO 11

LET H(I) = H(I) • 100 / H9

NEXT I

NEXT H
      610
                                                                                                                                                                                                                                                                                                                 5035 PRINT
5039 REMI INITIATE THE M N9 VERTICES
5040 LET N6 = 1
5050 LET N7 = 0
5060 LET N7 = 0
5060 LET N8 = 3
5070 PRINT "NUMBER OF PARAMETERS =";N8
5071 PRINT
5080 LET N9 = N8 + 1
5080 FOR I = 1 TO N9
5100 FOR I = 1 TO N9
5100 FOR I = 1 TO N9
5100 FOR J = 1 TO N8
5125 PRINT TWERTEX: ";I
5120 FOR J = 1 TO N8
5125 PRINT TWER (3);"PARAMETER # ";J;
5130 INPUT S(I,J)
5140 LET X(J) = S(I,J)
5150 NEXT J
                                                                                                                                                                                                                                                                                                                   5035
                                                                                                                                                                                                                                                                                                                                      PRINT
     630
640
650
      660
670
                          NEXT I
REM INITIATE HEADING AND PRINTOUT
       900
       905
                          PRINT
                        PRINT
PRINT
PRINT
PRINT
PRINT
TAB( 18);"C-13 NMR ANALYBIS OF"
PRINT TAB( 18);"MMA-Eta COPOLYMERS"
      906
      910
      920
                        PRINT TAB( 22); "SAMPLE : "; B$
PRINT TAB( 22): "SAMPLE : "; A$
      950
                          PRINT
PRINT "Y = MMA; Z = EA"

REM ANY COMPUTATIONS NEEDED CAN BE DONE HERE..
                                                                                                                                                                                                                                                                                                                                     NEXT J
GOSUB 3000
LET L(I) = I
                                                                                                                                                                                                                                                                                                                   5150
                                                                                                                                                                                                                                                                                                                   5160
      1000
                                                                                                                                                                                                                                                                                                                  5165
1001 REM ***
1999 GOTO 5000
2010 DATA "DBS", " B1 ", 20.5
2020 DATA "DBS", " B2 ", 18,2
2030 DATA "DBS", " B3 ", 16.5
2040 DATA "DBS", " A1 ", 50.2
2050 DATA "DBS", " A1 ", 50.2
2050 DATA "DBS", " A2 ", 51.6
2060 DATA "DBS", " A2 ", 51.6
2060 DATA "DBS", " C2 ", 176.7
2070 DATA "DBS ", "C2 ", 176.7
2090 DATA "DBS ", "C2 ", 176.7
2090 DATA "DBS ", "C2 ", 174.7
2110 DATA "DBS ", "C3 ", 174.7
2110 DATA "DBS ", "C5 ", 174.7
2110 DATA "DBS ", "C5 ", 174.7
2110 DATA "DBS ", "C4 ", 176.7
2000 REM THIS SUBROUTINE EVALUATES THE RESPONSE R9
3004 FOR K = 1 TO NB
3005 IF X(K) / 1.0 THEN 4700
3007 NEXT K
3010 LET G7 = X(1) + X(2)
3015 LET G9 = X(1) + X(2)
3015 LET G1 = 100 G7 GB GB / G9 / G9
3030 LET G(2) = 200 G7 * X(2) GB / G9 / G9
3030 LET G(3) = 100 * G7 * X(2) GB / G9 / G9
3050 LET G(5) = 100 * X(3)
3060 LET G(6) = 100 * X(3)
3071 LET G6 = 100 G7
X(3)
3080 LET G(8) = 2 GB * X(1) * X(2)
3080 LET G(9) = G6 * X(2) * X(2)
3081 LET G(9) = G6 * X(1) * X(2)
3081 LET G(9) = G6 * X(1) * X(2) * 2 * 6 6 * X(2) * X(3) * 4 * 6 6 * X(3)
3080 LET G(9) = G6 * X(1) * X(2) * 2 * 6 6 * X(1) * X(3) * 2 * 6 6 * X(3)
3080 LET G(9) = G6 * X(1) * X(1) * X(2) * 2 * 6 6 * X(2) * X(3) * 4 * 6 6 * X(3)
3080 LET G(9) = G6 * X(1) * X(1) * X(2) * 2 * 6 6 * X(2) * X(3) * 4 * 6 6 * X(3)
3080 LET G(9) = G6 * X(1) * X(1) * 2 * 6 6 * X(1) * X(3) * 2 * 6 6 * X(3)
3080 LET G(9) = G6 * X(1) * X(1) * 2 * 6 6 * X(1) * X(3) * 2 * 6 6 * X(3)
3080 LET G(9) = G6 * X(1) * X(1) * 2 * 6 6 * X(1) * X(3) * 2 * 6 6 * X(3)
   1001
1999
2000
                          REM ***
GOTO 5000
                                                                                                                                                                                                                                                                                                                   5170
                                                                                                                                                                                                                                                                                                                                          NEXT I
                                                                                                                                                                                                                                                                                                                  5170 NEXT ;
5171 PRINT
5172 PRINT
5180 PRINT "VX"; TAB( 5); "SX"; TAB( 10); "VTS RET/D";
5181 PRINT TAB( 26); "R VALUE"; TAB( 38); "PARAMETERS"
5182 PRINT
5190 FOR I = 1 TO N9
5200 PRINT 1; TAB( 4); NG; TAB( 25); 10 * R(I); TAB( 37);
5220 PRINT 1; TAB( 4); ";
5220 PRINT 1; TO N8
5230 PRINT 1; TAB( 4); ";
5240 PRINT 5(I, J); " ";
5240 PRINT 5(I, J); " ";
5240 PRINT 5(I, J); " ";
5245 PRINT
                                                                                                                                                                                                                                                                                                                                        PRINT

NEXT I

REM BEGIN SIMPLEX LOOP HERE

LET N6 = N6 + 1

GOSUB 4800

GOSUB 4900

GOSUB 4900
                                                                                                                                                                                                                                                                                                                    5260
                                                                                                                                                                                                                                                                                                                    527A
                                                                                                                                                                                                                                                                                                                  5280
5290
5300
                                                                                                                                                                                                                                                                                                                                          REM CALCULATE Y
                                                                                                                                                                                                                                                                                                                                    REM CRLCULATE Y
GOSUB 6300
FOR J = 1 TO N8
LET Y(J) = 2 ● D(J) - S(J9, J)
LET X(J) = Y(J)
LET S(10, J) = Y(J)
NEXT J
LET I = 9
CORUM 3000
                                                                                                                                                                                                                                                                                                                    5310
                                                                                                                                                                                                                                                                                                                   5340
                                                                                                                                                                                                                                                                                                                   5350
                                                                                                                                                                                                                                                                                                                  5380
                                                                                                                                                                                                                                                                                                                                         GOSUB 3000
LET J9 = 10
LET Z$ = "
                                                                                                                                                                                                                                                                                                                    5390
                                                                                                                                                                                                                                                                                                                                     GOSUB 4950
GOSUB 6300
GOSUB 6300
REM TEST FOR RESPONSE OF Y HERE
LET R1 = R(9)
IF R1 ( U(1) THEN 5500
IF R1 ) U(N9 - 1) THEN 5630
FOR J = 1 TO N8
LET S(J9, J) = Y(J)
NEXT J
LET R(J9) = R1
LET L(J9) = N7
GOTO 5260
REM CASE WHERE INTERVAL EXPANSION MAY BE NEEDED
LET M8 = 1
LET (C1 = 2
GOSUB 4750
LET 1 = 9
GOSUB 3000
                                                                                                                                                                                                                                                                                                                                         GOSUB 4950
                                                                                                                                                                                                                                                                                                                  5400
                                                                                                                                                                                                                                                                                                                  5404
5410
5420
5430
                       LET G(9) = G6 • X(1) • X(1) + 2 • G6 * X(1) * X(3) + 2 • G6 * X(3)
                       ● X(3)
                       • X(3)

LET G(10) = G6 • X(3) • X(3)

LET G(11) = 100 * X(3) * X(3) • X(3)

LET G1 = 0

FOR K = 1 TO V9

LET G1 = G1 + ABS (H(K) - G(K))
   3083
   3084
3090
                                                                                                                                                                                                                                                                                                                  5460
                                                                                                                                                                                                                                                                                                                 5470
5475
5480
5490
5500
   3100
   3110
                         LET G1 = G1 + ABS (H(K) - G(K))

NEXT K

LET R(I) = G1 / 110

RETURN

REM OPTIMIZATION OF MARKOVIAN PROBABILITIES (LL. 4000-6350)

REM FOR Y EXCEEDING LIMITS, RESET THE RESPONSE

LET R(I) = 10

RETURN

REM VECTOR ORDITION SUPPOLITIES Y-RESET X
   3120
3130
3300
                                                                                                                                                                                                                                                                                                                  5505
                                                                                                                                                                                                                                                                                                                  5510
5520
5530
   4000
   4700
4710
                                                                                                                                                                                                                                                                                                                 5530 LET I = 9
5540 GOSUB 3000
5550 LET R2 = R(9)
5560 IF R2 > = U(1) THEN 5450
5570 FOR J = 1 10 N6
5580 LET S(J9, J) = X(J)
5585 NEXT J
5586 ON MB GOTO 5587,5590,5592
5587 LET Z6 = "E"
                         REM VECTOR ADDITION SUBROUTINE Y=P+C*(P-X)

FOR J = 1 TO N8

LET X(J) = O(J) + C1 * (O(J) - S(J9, J))
    4750
   4760
4770
4780
4790
                          NEXT J
RETURN
                        NETURN
REM THIS SUBROUTINE SORTS THE RESPONSE R9
FOR I = 1 TO N9
LET M(1) = L(1)
LET U(1) = R(1)
NEXT I
FOR J = 1 TO N9 - 1
    4800
                                                                                                                                                                                                                                                                                                                 5587 LET Z* = _
5588 GOTO 5595
    4820
4825
```

```
GOSUB 4950
LET R (J9) = R2
LET L (J9) = R7
GOTO 5260
REM CASE WHERE INTERVAL CONTRACTION MAY BE NEEDED
IF R1 ( U (N9) THEN 5770
LET C1 = -0.5
GOSUB 4750
LET J = 9
GOSUB 3000
LET R3 = R (9)
IF R3 > U (N9) THEN 6000
FOR J = 1 TO NB
LET S (J9, J) = X (J)
NEXT J = "C"
GOSUB 4950
LET R (J9) = R (9)
LET L (J9) = R (9)
GOSUB 3000
LET R (J9) = R (9)
GOSUB 3000
LET R (J9) = R (9)
GOSUB 3000
LET R (J9) = R (9)
GOTO 5710
REM MARSIVE CONTRACTION APPLIED MERE
LET M8 = 2
LET C1 = .05
GOSUB 4750
LET L = .05
GOSUB 4750
LET L = .05
GOSUB 4750
LET L = .05
GOSUB 4750
LET C1 = .05
GOTO 5530
LET M = 3
LET C1 = .06
LET M = 3
LET C1 = .06
LET M = 7000
LET M = 7000
                                                                                                                                                                                                                                                                                                                                                                                  GOTO 6240
                                                                                                                                                                                                                                                                                                                                                                                GOTO 6240
RETURN
REM REPLACEMENT RESPONSE
FOR K = 1 TO N9
IF L(K) = M(N9) THEN 6340
NEXT K
LET J9 = K
RETURN
REM PRINTOUT OF FITTED RESULTS
PRINT
PRINT "NO" TOR( 11) *"I(ORS)": I
5660
5670
                                                                                                                                                                                                                                                                                                                                                                                   PRINT "NO"; TAB( 11);"I(OBS)"; TAB( 31);"I(CALC)"; TAB( 52);"DEV"
5691
5710
                                                                                                                                                                                                                                                                                                                                                                                   PRINT 1; TAB( 10);H(I); TAB( 31);G(I); TAB( 52);(H(I) ~ G(I))
NEXT I
PRINT
                                                                                                                                                                                                                                                                                                                                                           7150
7160
7170
5731
                                                                                                                                                                                                                                                                                                                                                                                    FOR K = 1 TO NB
PRINT " X(";K;") = ";X(K)
                                                                                                                                                                                                                                                                                                                                                            7180
                                                                                                                                                                                                                                                                                                                                                                                     PRINT " X(";K;") = ";X(K)
NEXT K
PRINT
REM CALCULATED COMPOSITION AND SEQUENCE DISTRIBUTION
REM ***
REM OTHER ADDITIONAL COMPUTATIONS MAY BE DONE HERE
REM END OF THE PROGRAM
PRINT "DO YOU HON'T TO DIEY MORE WITH THE DOTE 2".
                                                                                                                                                                                                                                                                                                                                                            7210
8000
                                                                                                                                                                                                                                                                                                                                                            8010
8028
9000
9910
                                                                                                                                                                                                                                                                                                                                                                                   PRINT
PRINT "DO YOU WANT TO PLAY MORE WITH THE DATA ?";
INPUT V$
IF V$ = "YES" THEN 5010
IF V$ = "NO" THEN 9970
GOTO 9930
PRINT "ANOTHER SAMPLE ?";
INPUT V$
INPUT V$
INPUT V$
IF V$ = "YES" THEN 9981
IF V$ = "NO" THEN 9998
GOTO 9971
RESTORE
                                                                                                                                                                                                                                                                                                                                                            9950
                                                                                                                                                                                                                                                                                                                                                            9960
                                                                                                                                                                                                                                                                                                                                                            9970
9971
9972
6104 GOSUB 4750
6105 GOTO 5530
6200 REM OPTIONAL TERMINATION
6210 IF N7 = 100 THEN 7000
6220 IF N7 - 10 * INT (N7 / 10) ) 0.1 THEN 6290
6230 PRINT "MORE? (VES/NO/RECAL)";
6230 IF V5 = "RECAL" THEN 5035
6260 IF V5 = "NO" THEN 7000
6270 IF V5 = "YES" THEN 6290
                                                                                                                                                                                                                                                                                                                                                            9973
                                                                                                                                                                                                                                                                                                                                                            9980
                                                                                                                                                                                                                                                                                                                                                                                        PRINT
                                                                                                                                                                                                                                                                                                                                                            9983
                                                                                                                                                                                                                                                                                                                                                                                      PRINT
                                                                                                                                                                                                                                                                                                                                                                                     GOTO 120
CALL EXIT
END
                                                                                                                                                                                                                                                                                                                                                            9984
```

Figure 3. Listing of the program FITCO, as applied to the EA-MMA copolymer case.

Table III. 13C Shifts of the Methyl Pentads of Polypropylene with Reaction Probabilities

¹³ C shift, ppm	assignment	first-order Markov	second-order Markov	bicatalytic sites		
21.8	mmmm	$P_{rm}(1-P_{mr})^3$	$\alpha^2 \gamma \delta$	$\omega(1-5\beta+5\beta^2)+(1-\omega)\sigma^4$		
21.6	mmmr	$(2P_{mr}P_{rm}(1-P_{mr})^2)$	$2\alpha\bar{\alpha}\gamma\delta$	$\omega(2\beta-6\beta^2)+2(1-\omega)\sigma^3(1-\sigma)$		
21.4	rmmr	$P_{rm}P_{mr}^{2}(1-P_{mr})$	$\bar{\alpha}^2 \gamma \delta$	$\omega\beta^2 + (1-\omega)\sigma^2(1-\sigma)^2$		
21.0	mmrr	$2P_{mr}P_{rm}(1-P_{mr})(1-P_{rm})$	$2\overline{\alpha\beta}\gamma\delta$	$\omega(2\beta-6\beta^2)+2(1-\omega)\sigma^2(1-\sigma)^2$		
20.8	mrmm + rmrr*	$2P_{mr}P_{rm}^{2}(1-P_{mr})$	$2\alpha\beta\gamma\delta$	$2\omega\beta^2+2(1-\omega)\sigma^3(1-\sigma)$		
		$\frac{+2P_{rm}P_{mr}^{2}(1-P_{rm})}{2P_{mr}^{2}P_{rm}^{2}}$	$+2\overline{\alpha\beta\gamma}\delta$	$+2\omega\beta^2+2(1-\omega)\sigma(1-\sigma)^3$		
20.6	mrmr	$2P_{mr}^2P_{rm}^2$	$2ar{lpha}etaar{\gamma}\delta$	$2\omega\beta^2+2(1-\omega)\sigma^2(1-\sigma)^2$		
20.3	rrrr	$P_{mr}(1-P_{rm})^3$	$\overline{lpha} \overline{oldsymbol{eta}} \delta^2$	$\omega\beta^2 + (1-\omega)(1-\sigma)^4$		
20.2	rrrm	$2P_{mr}P_{rm}(1-P_{rm})^2$	$2\overline{lphaeta}\deltaar{\delta}$	$2\omega\beta^2 + 2(1-\omega)\sigma(1-\sigma)^3$		
19.9	mrrm	$P_{mr}P_{rm}^2(1-P_{rm})$	$\frac{\alpha\beta\delta^2}{\alpha\beta+2\bar{\alpha}\delta+\gamma\delta}$	$\omega(\beta-3\beta^2)+(1-\omega)\sigma^2(1-\sigma)^2$		
	total	$P_{mr} + P_{rm}$	$\overline{\alpha\beta} + 2\bar{\alpha}\delta + \gamma\delta$	1.00		

Not resolved.




Figure 4. ¹³C NMR spectrum at 90 MHz of ethyl acrylate-methyl methacrylate copolymer.

intensities; thus, the smaller the R, the better is the fit.

The rest of the data in Table IV were taken from the heterogeneous catalysis work of Bukatov et al.²⁸ Sample B is the fraction soluble in boiling heptane, whereas sample C is the insoluble fraction. The soluble fraction B can be readily fitted to either first-order or second-order Markovian or bisite models, the latter giving only slight improvement. For the insoluble fraction C, either the second-order Markovian or the




Figure 5. ¹³C NMR spectrum at 90 MHz of the methyl region of polypropylene.

bisite model is preferred. The rather extreme values of the reaction probability parameters reflect the presence of steric propagation error, which generates only mmmr, mmrr, mrrm pentads.

As an example of how the FITCO program can be modified, the statements that need to be changed for polypropylene tacticity analysis are given in Figure 6. Otherwise the same program listing as in Figure 3 can be used.

3. Comonomer Sequence Determination. The use of computerized reaction probability models for sequence distributions has been previously reported for a number of copolymer

Table IV. Observed Intensities and Calculated Reaction Probabilities for the Methyl Pentads of Polypropylene

no.	pentad	sample A		sample B		sample C			
		obsd	calcd (M1)	obsd	calcd (M1)	calcd (BS)	obsd	calcd (M2)	calcd (BS)
1	mmmm	0.05	0.05	0.49	0.49	0.49	~0.91	0.91	0.91
2	mmmr	0.11	0.11	0.10	0.13	0.10	~0.05	0.04	0.04
3	rmmr	0.07	0.06	≤0.02	0.01	0.01	0	0.00	0.00
4	mmrr	0.11	0.12	0.11	0.10	0.11	~0.04	0.04	0.04
5	rrmr + mmrm	0.26	0.26	0.06	0.06	0.09	0	0.00	0.00
6	rmrm	0.12	0.14	≤0.01	0.00	0.03	0	0.00	0.00
7	rrrr	0.06	0.06	0.10	0.10	0.11	0	0.00	0.00
8	rrrm	0.15	0.13	0.06	0.08	0.07	0	0.00	0.00
9	mrrm	0.07	0.07	0.05	0.02	0.04	0	0.02	0.02
				First-Or	der Markov (M	1)			
	P_{mr}		0.54		0.118	-,		0.023	
	P_{rm}^{m}		0.50		0.295			0.995	
	R ^{''''}		0.0076		0.0122			0.0103	
				Second-C	order Markov (N	12)			
	P_{mmm}		0.449		0.882	,		0.980	
	P_{mrm}		0.462		0.302			0.010	
	P_{rmm}		0.478		0.877			0.992	
	P _{rrm}		0.513		0.295			0.991	
	R		0.0063		0.0122			0.0043	
				Bisit	e Model (BS)				
	$oldsymbol{eta}$				` /	0.058			0.018
	σ					0.249			0.428
	ω					0.669			0.998
	R					0.0096			0.0040

```
REM PROGRAM FITCO/PPTAC
REM DESIGNED FOR THE ANALYSIS OF METHYL PENTAD
REM OF POLYPROPYLENE AS OBSERVED VIA C-13 NMR
REM WRITTEN BY H. N. CHENG, HERCULES RESEARCH CENTER,
REM WILMINGTON, DELAWARE 19894, USA
                                   REM WILMINGTON, DELAWARE 19894, USA
REM PRINT "ANALYSIS OF METHYL PENTADS OF POLYPROPYLENE"
PRINT
PRINT
2 LET V9 = 9
3 READ V1$
3 PRINT I;"-"; TAB( 6);V1$;" AREA = ";
5 REM NORMALIZE THE INTENSITIES
5 LET H9 = 0
6 FOR I = 1 TO V9
6 FOR I = 1 TO V9
7 FOR I = 1 TO V
  500
550
560
570
                                                              PRINT TAB( 18); "ANALYSIS OF METHYL PENTAD OF PP"
                                                            PRINT

TAB( 18); "ANALYSIS OF METHYL PENTAD OF PP"

PRINT

REM ***

REM ANY COMPUTATIONS NEEDED CAN BE DONE HERE

REM ***

FOR K = 1 TO V9

IF H(K) ( ) 0 THEN 1006

LET H(K) = 0.00001

NEXT K

LET P1 = 2 ● H(1) / H(2)

LET E1 = P1 / (1 + P1)

LET P4 = 2 ● H(7) / H(8)

LET E2 = 1 / (1 + 2 ● P4)

LET P2 = H(4) / H(8)

LET E2 = 1 - P2 * E1

LET P3 = H(4) / H(8)

LET E3 = P3 ● (1 - E4)

PRINT "ESTIMATED 2ND ORDER MARKOV--Pmmm, Pmrm, Prmm"

PRINT TAB( 10);E1, "(";E2;")", "(";E3;")",E4

PRINT

GOTO 5000

REM THIS SUBROUTINE PROVIDES FORMAT FOR INPUT DATA

DATA "MMMM"

DATA "MMMR"

DATA "MMMR"
     1001
     1003
1004
1005
1006
                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                REM SIMPLEX SUBROUTINE
PRINT "INDICATE MODEL (1=MARKOV 1: 2=MARKOV 2: 3=BI-SITE)":
                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    5000
5035
     1010
                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           INPUT M9
PRINT "MARKOV 1--P(MM).P(RM): MARKOV 2--P(MMM).P(MRM),P(RMM),P(RRM
                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    5037
                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      PRINT "MARKOV 1--P(MM).P(RM): MARKOV 2--P(MMM).P(MRM),
)""
PRINT "BISITE MODEL---BETA,SIGMA,OMEGA"
REM INITIATE THE M N9 VERTICES
LET N6 = 1
LET N7 = 0
LET N8 = 3 THEN 5070
LET N8 = 3 THEN 5070
LET N8 = M9 * 2
REM CALCULATED COMPOSITION AND SEQUENCE DISTRIBUTION
REM ***
REM OTHER ADDITIONAL COMPUTATIONS MAY BE DONE HERE
REM ***
ON M9 GOTO 8660, 8690, 8740
PRINT "PRM =";X(2)
GOTO 8770
PRINT "PRM =";X(1)
PRINT "PRMM =";X(2)
PRINT "PRMM =";X(3)
PRINT "PRMM =";X(4)
BOTO 8770
PRINT "PRMM =";X(3)
PRINT "BETA = ";X(1)
PRINT "BETA = ";X(2)
PRINT "BETA = ";X(2)
PRINT "GMEGA = ";X(3)
PRINT "OMEGA = ";X(3)
     1040
1050
1060
                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 500440
500440
55560
500560
500660
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
606666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
606666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
606666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
606666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
606666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
60666
606
     1070
1080
1090
     1100
1110
1999
2000
  2010
2020
2030
                                                                                                                                                        "RMMR
                                                                        DATA
  2040
2050
2060
2070
                                                                        DATA
DATA
DATA
DATA
                                                                                                                                                        "MMRR
                                                                   DATA "RRRR"

DATA "RRRR"

DATA "RRRR"

DATA "MRRM"

ON M9 + 1 GOTO 9000, 3010, 3100, 3200

REM FIRST ORDER MARKOV

LET X(4) = X(2)

REM SECOND ORDER MARKOV

LET X(5) = 1 - X(1)

LET X(6) = 1 - X(2)

LET X(6) = 1 - X(2)

LET X(6) = 1 - X(4)

LET X(6) = 1 - X(4)

LET OP = 100 / (X(5) * X(6) + 2 * X(5) * X(4) + X(3) * X(4))

LET G(1) = 09 * X(1) * X(1) * X(3) * X(4) * X(5)

LET G(3) = 09 * X(3) * X(4) * X(5) * X(5)

LET G(4) = 09 * 2 * X(1) * X(3) * X(4) * X(5)

LET G(4) = 09 * 2 * X(1) * X(3) * X(4) * X(5)

LET G(4) = 09 * 2 * X(1) * X(6) * X(7) * X(8)

CETATORION CONTRACTOR OF THE PROPERTY OF THE PRO
                                                                                                                                                        "RRRR
     2080
2090
3009
       3010
                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    ALSO DELETE LINES 590-670;970;2100-2110;3020-3090;7190-7210 FROM THE FITCO/EAMMA CASE.
       3110
       3112
```

Figure 6. Statements needed in program FITCO to carry out model fitting in homopolymer tacticity.

systems, e.g., poly(ethylene-co-propylene), 12,13,17 poly(propylene-co-butylene),14 poly(ethylene-co-1-octene),15 and poly(propylene-co-1-octene). 15 It has been found that the computerized analytical approaches generally provide very precise values for the comonomer sequences even for copolymer spectra that contain overlapping resonances. 15,17

4. Chain Branching. The same approach has also been used for the determination of chain branching in low-density polyethylene.²⁹

USES AND LIMITATIONS

A major advantage of this computerized approach to polymer analysis is that all the NMR spectral lines are used for the calculations. Thus, experimental errors are minimized by this procedure. In addition, because the method checks for the intensities of all the resonance lines, if a wrong (or an uncertain) assignment is present, it often shows up when the analyses are carried out for a large number of samples. Thus, the model-fitting approach can assist in spectral interpretation.

In dealing with many types of polymers, this computerized approach has been found to be most valuable for the analysis of ¹H or ¹³C spectra with overlapping resonances. In these cases, although it is possible sometimes to devise computational schemes involving additions or subtractions of spectral intensities to obtain information of polymer composition or comonomer sequence distributions, in reality accurate numbers are often difficult to obtain, and if the overlap is severe, only estimates can be made. The use of the computerized reaction probability models thus can obviate these difficulties and provide precise results not otherwise available.

A final area of application is in the investigation of copolymerization mechanisms. Frequently one needs more information about a polymer system than simply the polymer composition and sequences. One may, for example, be interested in knowing the sensitivity of the reactivity ratios to changes in the initiators or the reaction conditions. One may be interested in testing the conformity of a polymerization process to a particular reaction probability model (e.g., Markovian or enantiomorphic). The tacticity of polypropylene is a good example. In these cases one can "squeeze" additional information from the NMR spectra with the present approach.

A limitation of the method is that the reaction probabilities are usually extracted from the NMR spectra of the final polymer product. As such, it only represents an average over the total polymerization process. If one is interested in the detailed description of the polymerization process as a function of conversion, one must repeat the analysis on polymer samples obtained at proper time intervals in the reaction vessel.

COMPARISON WITH SYNTHETIC APPROACH

In the synthetic approach proposed earlier, 10 one starts with a reaction probability model, obtains the expected sequence distributions, and then simulates a NMR spectrum. Comparisons can then be made visually of the experimental and the predicted spectra.

statistical model → sequence distribution → NMR shift behavior → NMR spectrum

The analytical approach reverses the process, starting with the spectrum and ending in the reaction probability parameters. The synthetic approach has the advantage that it is purely a computer experiment, with no polymer samples or spectrometer needed. Thus, different copolymers can be readily simulated. The analytical approach is much more rigorous but more time-consuming. It requires quantitative NMR spectrum with reliable spectral integrations and produces reaction probabilities that are most consistent with the spectral data at hand. The two approaches are actually complementary, and in many cases, the joint application of both approaches is desirable. In any unknown copolymer, for example, one can first use the synthetic approach to simulate a spectrum as an aid to spectral interpretation. Once the spectral lines are assigned, the analytical approach can be used routinely for the analysis of these copolymer samples.

EXPERIMENTAL SECTION

The program FITCO was written in BASIC for the Apple IIe. A listing of the program is shown in Figure 3. Alternate versions have also been made for Nicolet 1280 and PRIME 9950 computers. Interested readers may write to the author for copies of the Nicolet and the PRIME versions.

The ethyl acrylate-methyl methacrylate copolymers were either purchased (Monomer-Polymer, Trevose, PA) or made via free radical polymerization. These were dissolved as 15%(w/w) solutions in CDCl₃. The spectra were taken on a Nicolet NT360WB spectrometer at 90.55 MHz and ambient temperature. The free induction decays were stored in 16K memory addresses, zero-filled upon processing with a spectral window of 20 000 Hz. A pulse angle of 45° was used with 4-s delay. The T_1 values range from less than 0.2 s for the backbone carbons to approximately 0.6 s for the side chains. Combined with gated decoupling (thereby eliminating nuclear Overhauser effect), the conditions chosen should provide essentially quantitative spectra. Only carbonyl carbons give problems in quantitation. This has been taken care of in the program FITCO by separately normalizing all five carbonyl intensities to 100. (See Table II, footnote b.)

REFERENCES AND NOTES

- Bovey, F. A. High Resolution NMR of Macromolecules; Academic: New York, 1972.
- (2) Randall, J. C. Polymer Sequence Determination. Carbon-13 NMR
- Method; Academic: New York, 1977.
 (3) Pasika, W. M., Ed. Carbon-13 NMR in Polymer Science; ACS Symposium Series 103; American Chemical Society: Washington, DC,
- (4) Woodward, A. E., Bovey, F. A., Eds. Polymer Characterization by ESR and NMR; ACS Symposium Series 142; American Chemical Society: Washington, DC, 1980.
- (5) Ebdon, J. R. In Developments in Polymer Characterization; Dawkins, J. V., Ed.; Applied Science: New York, 1980; Vol. 2, pp 1-29.
- (6) Ernst, R. R. ACS Symp. Ser. 1982, No. 191, 47-61.
- Bax, A. Two-Dimensional Nuclear Magnetic Resonance in Liquids; D. Reidel: Dordrecht, Holland, 1984.
- Benn, R.; Gunter, H. Angew. Chem., Int. Ed. Engl. 1983, 22, 350-380. For example: (a) Bruch, M. D.; Bovey, F. A.; Cais, R. E. *Macromolecules* 1984, 17, 2547–2551. (b) Cheng, H. N.; Lee, G. H. *Polym*. Bull. (Berlin) 1985, 13, 549-556 and references cited therein.
- Cheng, H. N.; Bennett, M. A. Anal. Chem. 1984, 56, 2320-2327. Cheng, H. N. In Transition Metal Catalyzed Polymerizations. Alkenes
- and Dienes; Quirk, R. P., Ed.; Harwood Academic: New York, 1983; pp 617-638.
- Cheng, H. N. Anal. Chem. 1982, 54, 1828-1833.
- (13) Carman, C. J.; Harrington, R. A.; Wilkes, C. E. Macromolecules 1977, 10, 536-544.
- Cheng, H. N. J. Polym. Sci., Polym. Phys. Ed. 1983, 21, 573-581.
- (15) Cheng, H. N. Polym. Commun. 1984, 25, 99-105.
 (16) Cheng, H. N. Polym. Bull. (Berlin) 1985, 14, 347-354
- Cheng, H. N.; Lee, G. H. Polym. Bull. (Berlin) 1984, 12, 463-470.
- Bovey, F. A. High Resolution NMR of Macromolecules; Academic: New York, 1972.
- (19) Koenig, J. L. Chemical Microstructure of Polymer Chains, Wiley-Interscience: New York, 1980.
- (20) Sheldon, R.; Fueno, T.; Furukawa, J. J. Polym. Sci., Polym. Phys. Ed. 1969, 7, 763-773 and references therein.
- (21) Zambelli, A.; Locatelli, P.; Provasoli, A.; Ferro, D. R. Macromolecules **1980**, 13, 267-270.
- Schilling, F. C.; Tonelli, A. E. Macromolecules 1980, 13, 270-275.
- Cheng, H. N.; Lee, G. H. Macromolecules, in press. Zhu, S.-N.; Yang, X.-Z.; Chujo, R. Polym. J. (Tokyo) 1983, 15, 859-868.
- (25) Inoue, Y.; Itabashi, Y.; Chujo, R.; Doi, Y. Polymer 1984, 25, 1640-1644
- Pavan, A.; Provasoli, A.; Moraglio, G.; Zambelli, A. *Makromol. Chem.* 1977, 178, 1099–1109. (26)
- Suter, U. W.; Neuenschwander, P. Macromolecules 1981, 14, 528-532. Bukatov, G. D.; Zakharov, V. A.; Yermakov, Y. I.; Zambelli, A.
- Makromol. Chem. 1978, 179, 2093-2096.
- (29) Cheng, H. N. Polym. Bull. (Berlin) 1986, 16, 445-452.