


**...you asked for
reproducible,
straight line,
separations
...we gave you**



CASCADE ELECTRODE CURTAIN ELECTROPHORESIS

...widely used for separations of:

- Viruses and Nucleic Acids •
- Hormones • Immunizing Agents •
- Antibiotics • Inks, Stains and Dyes •
- Bacterial Extracts • Amino Acids •
- Enzymes • Cell Fractions • Pla-
- cental Fractions • Spinal Fluid •
- Cell and Tissue Homogenate
- Extracts • Blood Proteins, Lipopro-
- teins and Mucopolysaccharides •
- Inorganics (e.g., Noble Metals,
- Rare Earths and Nuclear Fission
- Products) • Botanical Products.

— FEATURES —

- Greater speed, accuracy and reliability • Minimum heat generation, no external cooling required • Low ripple, constant current power supply • Uniform spreading of buffer • Entire electro-conductive surface is available for migration of fractions on the paper curtain.

CASCADE ELECTRODE ELECTROPHORESIS

...gives straight line separations when used to perform: • Fraction Collections—isolates material quantitatively • Continuous Flow Electrophoresis • Preparative Strip Chromatography • Star Electrophoresis • Electrochromatography.

— REFERENCES —

Hub. Peeters and W. de Keersgieter, **ELECTROCHROMATOGRAPHY OF THE LIPIDS OF NORMAL SERUM**. *Prot. of the Biol. Fluids*, 1961, 225-232

K. Maertens and H. Peeters, **ELECTROCHROMATOGRAPHY OF HUMAN OCULAR LENS TISSUE**. *Prot. of the Biol. Fluids*, 1961, 253-262

Sobtko, H. & Stewart, C. P. **ADVANCES IN CLINICAL CHEMISTRY**, VOL. 2, 1959, 91-134

AT YOUR LABORATORY SUPPLY DEALER
another quality product of


**Scientific
Industries INC.**



Dept. AC664, 220-05 97th AVENUE
QUEENS VILLAGE 29, NEW YORK

Circle No. 156 on Readers' Service Card

SMALL BOMBS and PRESSURE VESSELS



Sizes: 22, 45, 71, 128, 200, 310, 510 ml, and larger.
Materials: T316 stainless, Monel and other alloys.
Working pressures: up to 5000 psig. and higher.
Maximum temperatures: 350°, 500° C. and higher.
Valves, electrodes and other fittings can be provided.
Let PARR help you select a combination suitable for your particular requirements.

PARR INSTRUMENT COMPANY

211 Fifty-Third St.

Moline, Illinois

Circle No. 88 on Readers' Service Card

1963 DIRECTORY OF GRADUATE RESEARCH

The newest edition of this unique directory is the sixth to be prepared by the ACS Committee on Professional Training. It covers the 1961-62 and 1962-63 academic years and provides a useful reference to:

- degrees available
- fields of interest and publications of 4,152 faculty members

in 297 departments or divisions of chemistry, biochemistry, and chemical engineering in United States universities offering the Ph.D. degree.

Under each department heading, degrees offered and fields of specialization appear first. Then faculty members are listed alphabetically with an up-to-date record on their education . . . general fields of major research interest . . . subjects of current research . . . publications during the past two years.

The table of contents lists universities under the three main groups. There is another index by faculty names. A summary table shows for each graduate department of chemistry the number on the faculty, number of Ph.D.'s in each department, graduate enrollment, and Ph.D. degrees granted in 1961-62 and 1962-63. It offers a concise comparison by size.

If you counsel students or seek an advanced degree yourself, or if you are interested in knowing the kind of research done in certain academic centers for whatever purpose, then this book will answer your questions and save you time.

655 pages.

Paper bound.

Price: \$4.00

Order from:

Special Issues Sales American Chemical Society
1155 Sixteenth Street, N.W., Washington, D. C. 20036