



Nine consecutive injections of C14, C15, and C16 hydrocarbons.

NO LOSS IS YOUR GAIN.

Nine 1 μ L injections from only 10 μ L of sample.

Imagine being able to perform multiple injections with such a minute amount of sample! Varian does. With our new 8100 AutoSampler™, we have redefined low loss autosampling into virtually no loss.

Imagine the possibility of no sample carryover! None whatsoever. The 8100 AutoSampler is unique because of the dynamic purge design of the wash system. You inject the whole sample, and nothing but the sample.

The impossible becomes possible because of the superior capabilities of the 8100 achieved with the "sandwich" technique — a technique that not only improves precision and accuracy, but also ensures that what is measured is injected.

With an 8100 installed on your Varian Gas Chromatograph, you'll easily see the benefits to your lab in increased productivity. It's more than just a sample saver. It's easy to use because it's completely controlled via the keyboard on the GC.

AutoSampler methods are stored and battery protected, so for unattended operation, just load it and go.

Let us show you how the "no loss" 8100 AutoSampler can be your gain. For more information, call **1-800-926-3000**. In Canada, call **416-457-4130**.

CIRCLE 140 ON READER SERVICE CARD

GC WITH A FUTURE

varian 