

During our visits to the field laboratories and in our communications with them, we frequently hear, "We don't have time to spend on all of this quality assurance work. It cuts sharply into our production". The answer is, of course, "we do not have time *not* to perform the quality assurance checks". Of what value is it to turn out reams and reams of data if we do not have the assurance, other than a gut feeling—assurance with facts and figures that the data obtained are accurate and reliable?

In addition to the above, another safeguard FDA applies to its use of analytical methods is the requirement that if a product is found by analysis to be in violation of the FD&C Act, a second or "check" analysis has to be performed. The check analysis has to be run by a second analyst. The check

analyst starts from scratch making his own composite, if one is required, checking his own standards, solutions, etc. Only if the results from the two analyses are in reasonable agreement is regulatory action considered against the product. Any method used which is not official must be validated in the hands of the analyst by recovery studies or other appropriate studies before the results by the method are acceptable. We take into account the fact that frequently such a quantitative method is not absolutely definitive, and we must follow the quantitative analysis with a confirmative step such as mass spectrometry or thin-layer chromatography. We must be as certain as possible that the ingredient we are quantitating is in fact the ingredient we think it is.

So goes the battle in FDA's regulatory laboratories.

New Approaches to FDA Analytical Problems[†]

THOMAS CAIRNS* and ROBERT A. JACOBSON

U.S. Food and Drug Administration, Los Angeles, California 90015

Received January 14, 1977

The formulation of a molecular weight listing of pesticides and industrial chemicals is described. Application to residue analysis by mass spectrometry is reported. A novel approach to quantitation of polychlorinated biphenyls is presented.

INTRODUCTION

During normal surveillance analysis for pesticide residues, compounds are frequently encountered which cannot be readily identified by available gas chromatographic data. In the past, identification was attempted using different detectors and stationary phases in the gas chromatographic separation in the hope that the combined retention data acquired might correspond to an already characterized pesticide or industrial chemical. Such cross correlations of retention data are time-consuming and rely heavily on the recreation of the standard conditions under which the original data for comparison was formulated. Furthermore, the compilation of a substantial retention data base for this type of approach to identification is a long-range undertaking to achieve absolute completion.

BACKGROUND

More recently, with the acquisition of combined gas chromatography-mass spectrometry (GCMS) systems in seven FDA field laboratories, the structural elucidation of detected unidentified compounds by conventional GC has been referred directly to GCMS techniques. The essence of this approach to structural elucidation, especially with low-resolution instruments, is knowledge of the molecular weight of the compound under investigation. In particular, chemical ionization (CI) techniques often provide conclusive evidence for the recognition of the so-called quasi-molecular ion or protonated molecular ion. This single piece of information alone can narrow down the possibilities provided a comprehensive listing of pesticides and industrial chemicals is available, organized according to molecular weight. No commercially available compilation¹ exists which contains within its database all the compounds pertinent to FDA analyses.

MOLECULAR WEIGHT LISTING OF PESTICIDES AND INDUSTRIAL CHEMICALS

The need for a listing containing all the possible compounds

DATE - 8-28-76	COMPOUND NAME	FORMULA	PAGE 1
51 0265	ACRYLONITRILE	C ₃ H _{3.5} N	REN
56 0262	ACROLEIN	C ₃ H ₄ O	ACV
58 0419	ALLYL ALCOHOL	C ₃ H ₆ O	ALP
62 0368	ETHYLENE GLYCOL	C ₂ H ₆ O ₂	
64 0007	SODIUM CYANAMIDE	C ₂ H ₂ N ₂ NA	0002
73 0891	SEC-BUTYLAMINE	C ₄ H ₁₁ N	0482
76 0617	B-HYDROXYETHYLHYDRAZINE	C ₂ H ₆ N ₂ O	0078V2
77 0277	FLUORACETAMIDE	C ₂ H ₄ F N O	FLP
84 0416	ANITROLE	C ₂ H ₄ NA	0004AMT
86 0268	2-BUTYNE-1,4-DIOL	C ₄ H ₈ O ₂	BYV
88 0524	ETHYL ACETATE	C ₄ H ₈ O ₂	ETP
92 0829	1-CHLORO-2,2-DIETHOXYPROPANE	C ₅ H ₁₀ Cl O ₂	EP1
93 0622	CHLORACETIC ACID	C ₂ H ₂ Cl O ₂	CPA
94 0531	HYTROL 100	C ₅ H ₈ N ₂	0474HYN
97 0640	1,2-DICHLOROETHANE	C ₂ H ₄ Cl ₂	
99 0143	ALLYL ISOTHIOCYANATE	C ₄ H ₅ N S	ALC
99 0916	SODIUM MONOFLUOROPHOSPHATE	C ₂ H ₂ F O ₄ NA	0275HF
102 0252	ETHYLENE DITHIOUR	C ₂ H ₆ N ₂ S	0000H1
105 0822	4-CHLOROACRYLIC ACID	C ₃ H ₃ Cl O ₂	CH
106 0419	BENZALDEHYDE	C ₇ H ₆ O	BEV
107 0908	3-CHLOROPROPIONIC ACID	C ₃ H ₅ Cl O ₂	CP
109 0690	1,2-DICHLOROPROPENE (D-)	C ₃ H ₄ Cl ₂	0005D1P
110 0268	ORTHOHOL, O-DIHYDROBENZENE	C ₆ H ₆ O ₂	COL
111 0085	ARINOMETHYLPHOSPHONIC ACID	C ₆ H ₅ N O ₃ P	0085
111 0614	N-VINYL-2-PROPIONOLONE	C ₅ H ₈ N O	1204
111 0647	1,2-DICHLOROPROPANE (D-)	C ₃ H ₆ Cl ₂	0005D1P
112 0080	CHLOROBENZENE	C ₆ H ₅ Cl	0245
112 0272	MALEIC HYDRAZIDE	C ₄ H ₄ N ₂ O ₂	0005MHV
114 0099	ARMINOL SULFONATE	HS N ₂ O ₃ S	0000
120 0067	FURAN	C ₄ H ₄ O	0181F0P
122 0047	DICHLOROBUTENES	C ₄ H ₆ Cl ₂	0181F0P
128 0626	NAPHTHALENE	C ₁₀ H ₈	0181F0P
128 0680	UREAM, METHAN-SODIUM	C ₂ H ₄ N ₂ SO ₂ NA	0150MEU
128 0901	4-CHLOROPYRIDINE-N-OXIDE	C ₅ H ₄ Cl N O	CP1
129 1156	ISOBUTYL ALCOHOL	C ₄ H ₁₀ O	1021
131 0508	1,4,1-TRICHLOROETHANE	C ₂ H ₃ Cl ₃	0193TR1
135 0040	2-METHOXYNITROBENZOLE, METHYL OF BENZYL	C ₇ H ₇ N ₂	0411BE1
135 0524	EPICLOMETHYLIN	C ₅ H ₈ Br O	1205
137 0662	CHLORACILIC ACID	C ₂ H ₂ O ₂ AS	0400CAC
139 0455	METHYL ARSONIC ACID	C ₂ H ₅ O ₂ AS	0420MAA
141 0012	MUNITION	C ₂ H ₅ N O ₂ P S	0420MAA
141 0185	1-BROMO-2-CHLOROETHANE	C ₂ H ₄ Br Cl	ETN
141 0588	CHLORON	C ₂ H ₄ Cl ₂ O ₂	0000H1
141 0952	CHLORON, BIS-(2-CHLOROETHYL)ETHER	C ₄ H ₈ Cl ₂ O	0180CH
142 0185	4-CHLORO-O-CRESOL, METABOLITE OF MCP	C ₇ H ₇ Cl O	0562
142 0185	ETHIDE, 1,1-DICHLORO-1-NITROETHANE	C ₂ H ₃ Cl ₂ N O ₂	0005EDT
143 0742	ETHYLPHON	C ₂ H ₅ Cl O ₂ P	0410ETH
145 0091	FRYCHLOR HYDROLYSIS PRODUCT	C ₁₀ H ₁₁ N	0459
145 0690	M-DICHLOROBENZENE	C ₆ H ₄ Cl ₂	0244
145 0690	O-DICHLOROBENZENE	C ₆ H ₄ Cl ₂	0190DBE
145 0690	P-DICHLOROBENZENE	C ₆ H ₄ Cl ₂	0191PBD
146 0044	ACETONE, PHENYLHYDRAZONE	C ₉ H ₁₀ N ₂	ACZ
147 0610	2-ETHYL-1,3-HEXANEDITHIOL	C ₈ H ₁₈ O ₂	EXU
147 0610	CARYONE, F-METHANE, S-DIEN-2-ONE	C ₁₀ H ₁₄ O	CCP
148 0088	ANITROLE, 1-METHOXY-4-PROPENYLBENZENE	C ₁₀ H ₁₂ O	ANE
149 0607	CORH	C ₆ H ₁₂ Cl N O	0189COR

Figure 1. A typical page from the Molecular Weight Listing of Pesticides and Industrial Chemicals.

that might be encountered in regulatory analyses is long overdue. To meet this demand, such a database has been initiated and at the moment the listing containing 1650 compounds has been computer sorted according to molecular weight (Figure 1). The identification code refers to the




Figure 2. GCMS residue analysis of cherries.

"Primary Standard Data Sheet" provided with each standard supplied to FDA personnel. The three letter code appearing after the numerical code is from the "Nanogen Index".² Both these codes are references which provide additional details on structure, toxicity, usage, etc. Although the major anticipated use of this list involves mass weight to the nearest integral whole number, exact molecular weights to the fourth decimal

place have been deliberately incorporated to aid in high resolution measurements and mass defect determinations. As each of these compounds or their metabolites are recorded by an FDA mass spectrometrists, the resulting mass spectral data (EI, CI, FI, or FD) will be numerically included (ten strongest peaks and intensities) in the listing and concurrently printed up as a plotted normalized spectrum for a master compilation of reliable mass spectral data. It is planned to develop further indices to this master compilation by name and molecular formula. The unique feature of this database of mass spectral information which differentiates it from any other is that the pesticides and industrial chemicals to be encountered have been previously identified and sorted by molecular weight without waiting for the characteristic mass spectrum.

EXAMPLES

Case 1. Residue Analysis of Cherries. To illustrate the potential of such a generated database, a recent example of an unidentified peak or eluting compound from the residue analysis of fresh Bing cherries was quickly solved by reference to the molecular weight listing of pesticides and industrial chemicals. The total ion chromatogram (TIC) (Figure 2) indicated a number of peaks in the sample. However, by the simple computer technique of searching each mass scan for only m/e 149, an ion highly characteristic of phthalates, it was clear that peaks a, b, c, and d belonged to this class of compound. The resulting mass spectrum of the unknown revealed that it had a molecular weight of 367. Positive identification of the molecular weight was observed by the presence of an adjunct ion at m/e 396 corresponding to $(M + 29)^+$ resulting from an ion-molecule collision with $C_2H_5^+$ formed by high-pressure methane collisions. In addition, the presence of an ion at m/e 370 with approximately 30% of the intensity of the MH^+ ion at m/e 368 was strong indication that the molecule possessed one chlorine atom. A quick visual search of the molecular weight listing indicated that the

MOLECULAR WEIGHT LISTING OF PESTICIDES AND INDUSTRIAL CHEMICALS			
DATE - 8-16-76		PAGE 17	
NO.	COMPOUND NAME	FORMULA	ICENT
151 8550	1-HYDROXYCHLOROBENZENE	C18 H6 CL6 O	8405
152 8447	O,P-DITHIOPHOSPHATE	C14 H8 CL5	8201000
153 8447	P,P'-DDT	C14 H9 CL5	820000P
154 8840	TETRAPHENYL	C12 H6 CL4 O2 S	812870
155 8840	BITHIOPHOSPHATE	C12 H6 CL4 O2 S	812870
156 8332	2,4,5-T-BUTYL ETHYL ESTER	C14 H17 CL2 O4	821107
157 8427	CARBOPHOSPHATE-1,1,1,1-TETRAPHENYL ETH	C16 H22 N2 O7	8552
158 8427	METHYL DINITRO-OCTYL-PHENYL CARBONATES	C16 H22 N2 O7	8552
159 8517	FLUOROPHOSPHATE	C12 H17 CL F3 N2 O4	8517FIN
160 8517	ACETOPHOSPHATE	C14 H14 CL N2 O2 P S	8517P
161 8517	ACETOPHOSPHATE	C12 H14 O4	1194
162 8444	HEXACHLORO BIPHENYLS (PCBS)	C12 H4 CL6	811
163 8444	CARBOPHOSPHATE SULFOXIDE	C11 H16 CL O3 P S2	811
164 8444	BETA (N-1,1,1,1-CHLOROPHENYL)PHOSPHATE	C12 H14 CL2 O4 P	8459CHM
165 8444	ALPHA (N-1,1,1,1-CHLOROPHENYL)PHOSPHATE	C12 H14 CL2 O4 P	8459CHM
166 8444	CARBOPHOSPHATE OXYGEN ANALOG SULFONE	C11 H16 CL O5 P S2	8271
167 8444	TRITHIOPHOSPHATE GLYCOL DIBENZOTRIS	C20 H22 O6	1146
168 8444	ETHYL ANTI-OXIDANT 706	C22 H30 O2 S	1804
169 8444	CHLOROPHOSPHATE	C11 H15 CL2 O2 P S2	811
170 8444	URACIL MONOPHOSPHATE-TCA	C11 H12 CL4 N2 O6	8099H0T
171 8444	AVENGE	C18 H20 N2 O4 S	8556
172 8444	2,4-DIB. ISOOCTYL ESTER	C18 H26 CL2 O2	8197
173 8444	TETRAETHYL CITRATE	C18 H22 O7	1194
174 8444	ISOCIN	C12 H8 CL6	80941SN
175 8444	AMEN. ALERIN	C12 H8 CL6	80941SN
176 8444	CUMAPHOS	C14 H16 CL O5 P S	801500P
177 8444	2,4-D. BEP ESTER	C17 H24 CL2 O4	8226
178 8444	2-ETHYLHEXYL DIPHENYL PHOSPHATE	C20 H27 O4 P	1228
179 8444	BENZONATE	C18 H18 CL N O5	8197
180 8444	BENZOPHOS	C8 H8 BP CL2 O2 P S	84298SP
181 8444	GARDON. TETRACHLOROPHOSPHATE	C10 H9 CL4 O4 P	81887EP
182 8444	EPHON	C11 H9 CL5 O2	85288EP
183 8444	TETRACHLOROETHYLPHENOL A	C15 H12 CL4 O2	1042
184 8444	O,O-DIETHYL O-NAPHTHYLAMINO PHOSPHOROTHIO	C16 H16 N O5 P S	811
185 8444	SUFFIN	C18 H17 CL2 N O3	857980H
186 8444	CHLOROCYCLOHEXANE	C7 H12 CL2	811
187 8444	HEPTACHLOROPHTHALENE	C10 H4 CL7	811
188 8444	TETRACHLORO TERPHENYLS	C18 H10 CL4	811
189 8444	ANISURON	C17 H16 CL2 N2 O2	811
190 8444	2,4,5-T. ISOOCTYL ESTER	C18 H21 CL3 O2	822700
191 8444	2,4,5-T. 2-ETHYL HEVYL ESTER	C18 H21 CL3 O2	822700
192 8444	2,4,5-T. 2-ETHYL HEVYL ESTER	C20 H20 O6	1182
193 8444	PHOSALONE	C15 H15 CL N O4 P S2	8199PHE
194 8444	OC CHA C I		
195 8444	184 182 155 182 186 170 165 137 154 152		
196 8444	48 28 28 26 25 18 16 9 8		
197 8444	O,P'-HEPTHANE	C14 H9 CL5 O	8442
198 8444	O,P'-HEPTHANE	C14 H9 CL5 O	8442
199 8444	QUINOSIN C	C6 H11 CL O2 H6	811
200 8444	TPH. FENTIN HYDROXIDE	C15 H16 O SN	8156FEL
201 8444	2,4,5-T PROP. GLY. BU. ETHER ESTERS	C15 H19 CL3 O4	82107CP
202 8444	TRICRESYL PHOSPHATE	C21 H21 O4 P	1226
203 8444	GLYCOLIN	C22 H44 N2 O2	80956LO
204 8444	4,5-DIBROMOSALICYLANILIDE	C17 H9 BR2 N O2	1207
205 8444	HEPTACHLOR	C10 H5 CL7	8082HEP
206 8444	DINOSULFON	C16 H22 N2 O6 S	1201
207 8444	TOXINIL	C7 H2 N2 O	8517

Figure 3. Identification of Phosalone and primary standard data sheet referred to from listing.

FDA No. 199	
PRIMARY STANDARD	
DATA SHEET	
FOOD AND DRUG ADMINISTRATION	
Pesticides Reference Standards Section, Pesticides Branch, Division of Food Standards and Additives, Bureau of Science, Washington, D.C. 20204	
FDA No. 199	
Name: phosalone (USASI name); Zalone*	
Systematic Name: O,O-Diethyl S-(6-chloro-2-oxobenzoxazolin-3-yl) methyl phosphorodithioate	
Empirical Formula: C ₁₅ H ₁₅ ClN ₂ O ₄ P ₂ S ₂	M.W. 367.82 States: Solid
Structure:	
Stability: Stable	
Toxicity: Acute oral LD ₅₀ (rats) 135mg/kg. Absorbed through the skin.	
Purity: 99.12 ± 0.52	
Differential Scanning Calorimetry: 99.12 ± 0.17 mole percent	
Gas Chromatography, thermionic flame:	
One peak, rel. ret. time (parathion) 5.50	
Response: 26ng equivalent to 2ng parathion	
Conditions: 6' 10% DC 200 on 80/100 Gas Chrom Q, column 200', 60 ml/min N ₂ , 0.2 x 10 ⁻⁸ amps	
Melting Range, capillary, 1°/min: 45.6 - 47.7° C.	
July 1967	




Figure 4. GCMS residue analysis of fresh carp.




Figure 5. Identification of DDMS and plotted normalized mass spectrum for inclusion in the FDA Mass Spectral Data Compilation of Pesticides and Industrial Chemicals.




Figure 6. GCMS analysis of Aroclor 1254 and subsequent limited mass scan searches for chemical profiling.




Figure 7. GCMS analysis of salmon extract illustrating co-eluting compounds at scan 212.

compound was Phosalone (Figure 3). Unambiguous identification was provided by recording a standard reference sample under exactly similar conditions as the sample. Both the retention time and mass spectrum agreed. The unknown




Figure 8. GCMS analysis of salmon extract and limited mass scan searches to profile PCB's and DDE.




Figure 9. GCMS-CI-SIM analysis of salmon extract.

had been quickly solved. A point worthy of mention is that a solution was arrived at by a low-resolution instrument. Normally, such a complex molecule containing N, O, S, P,




Figure 10. GCMS-CI-SIM analysis of Aroclor 1254.

and CI would have required high-resolution measurements to suggest plausible structures.

Case 2. Residue Analysis of Fresh Carp. Sometimes the mass spectrum of a compound under investigation has not yet been added to the database although the molecular weight, molecular formula, etc., have already been added. Such information can provide clues to the possible identity of the unknown. To illustrate this additional facet of the uncompleted database, an unidentified compound was encountered in the analysis of fresh carp. A number of peaks had already been identified by existing GC using electron capture detector (EC) systems. In this particular case the unknown was strongly suspected of being a metabolite both by its retention time and close similarity in mass spectral characteristics to the identified metabolites of DDT found present, DDMU and DDD (Figure 4). By inference, the molecular weight was calculated to be 284 since the prevailing fragmentation of this type of structural entity was to lose a chlorophenyl group, i.e., loss of 111 amu. A search of the listing (Figure 5) revealed that the possible metabolite could be DDMS. Inspection of an available standard of DDMS proved the case and the spectrum was then added to the list and compilation for future reference.

In summary, the construction of a master listing of pesticides and industrial chemicals used today together with their respective mass spectra can provide a powerful tool to aid in the structural elucidation of unknowns found in regulatory samples. If the compound is not on the list then at least a number of possibilities have been instantly removed from consideration and further detailed study will be necessary to solve the identity of the compound. This screening process provides either quick solutions or difficult problems requiring extended study enabling the analyst to make decisions on time and effort to be expended.

CHEMICAL PROFILING OF POLYCHLORINATED BIPHENYLS

Quantitative analysis of polychlorinated biphenyls (PCB's) has been a difficult problem since their detection in the environment by Jensen³ in 1966, primarily because of the large number of compounds involved. A novel approach recently experimented with PCB residue samples has been the chemical profiling of the PCB content in the sample to ascertain the

most suitable available standard to be used for subsequent quantitation even in the presence of interfering substances. The essence of the success of this approach was based on the observed fact that PCB's on CI produce almost entirely the molecular ion cluster. Since only molecular ions are formed, limited mass scans (LMS) can be conducted using the data system to profile the PCB content by specifying those ions corresponding to monochloro, dichloro, trichloro, etc. The resulting LMS searches were found to be highly characteristic in differentiating the various Aroclor standards available (Figure 6) and could be effectively used as chemical profiles to choose which Aroclor standard was most suitable for quantitation of a residue sample. To illustrate the real world problems encountered in the analysis of PCB's in fish samples (Figure 7), the mass spectrum represented by scan 212 in the TIC chromatogram indicated that this particular peak contained two co-eluting compounds—DDE and a pentachlorobiphenyl. Such interferences make selection of a suitable standard for quantitation difficult and impose further chemical separation and clean-up procedures. However, by computer techniques on this TIC chromatogram, LMS searches were able to pull out selectively both the PCB and DDE profiles (Figure 8). The observed profile for the salmon extract was similar to that previously observed for Aroclor 1254. Certain discrepancies existed in the profiles which could only have resulted from weathering and/or metabolism.⁴ However, the overall profile closely resembled the profile observed for Aroclor 1254. To achieve quantitation without further chemical treatment of the sample the extract was rerun using single ion monitoring (SIM) techniques (Figure 9) and the areas under each ion value were calculated and compared with those observed for a standard Aroclor 1254 sample (Figure 10). The calculated value of 14 ppm agreed with the value using conventional GC with EC.

CONCLUSIONS

With the experience and positive reactions gained from the molecular weight listing project it is anticipated that additional compounds and metabolites will be added to the list as they are reported in the literature. A deliberate attempt to record all the mass spectral data is currently underway as a joint effort within FDA laboratories. At the same time, available mass spectral data are being forwarded to EPA for inclusion in the Cyphernetics MSSS file system.¹

In the case of PCB quantitation, plans are already underway to refine the chemical profiling techniques and quantitate on the basis of individual isomeric groups (mono, di, tri, etc.) rather than related to an Aroclor standard mixture.

ACKNOWLEDGMENT

The authors wish to express particular thanks and appreciation to Jerry Froberg for helpful discussions and suggestions.

LITERATURE CITED

- † Presented before the Division of Chemical Information, 172nd National Meeting of the American Chemical Society, San Francisco, Calif., Aug. 31, 1976.
- (1) R. S. Heller, G. W. A. Milne, R. J. Feldmann, and S. R. Heller, "An International Mass Spectral Search System (MSSS). V. A Status Report", *J. Chem. Inf. Comput. Sci.*, **16**, 176 (1976).
- (2) K. Packer, "Nanogen Index. A Dictionary of Pesticides and Chemical Pollutants", Nanogens International, Freedom, Calif. 95019, 1975.
- (3) S. Jensen, "Report of a New Chemical Hazard", *New Sci.* **32**, 612 (1966).
- (4) L. D. Sawyer, "Collaborative Study of the Recovery and Gas Chromatographic Quantitation of Polychlorinated Biphenyls in Chicken Fat and PCB-DDT Combination in Fish", *J. Assoc. Off. Anal. Chem.*, **56**, 1015 (1973).